








Tabeeb Angiography Catheter						
Tip Configuraion	Size (Fr)	Length (cm)	Accept Guidewire OD (in)	Side holes	Shape	Applicable Vessels
JL 3.5	5	100	0/038	0		Left Coronary
JL 3.5	6	100	0/038	0		
JL 4.0	5	100	0/038	0		
JL 4.0	6	100	0/038	0		
JL 5.0	5	100	0/038	0		
JL 5.0	6	100	0/038	0		
JR 3.5	5	100	0/038	0		Right Coronary
JR 3.5	6	100	0/038	0		
JR 4.0	5	100	0/038	0		
JR 4.0	6	100	0/038	0		
JR 5.0	5	100	0/038	0		
JR 5.0	6	100	0/038	0		
TIG 3.5	5	100	0/038	0		Left and Right Coronary
TIG 3.5	6	100	0/038	0		
TIG 4.0	5	100	0/038	0		
TIG 4.0	6	100	0/038	0		
MPA	5	100	0/038	0		Multi - Purpose
MPA	6	100	0/038	0		
MPA1	5	100	0/038	0		
MPA1	6	100	0/038	0		
MPA2	5	100	0/038	0		
MPA2	6	100	0/038	0		
AL1	6	100	0/038	0		Left and Right Coronary
AL2	6	100	0/038	0		
AR1	6	100	0/038	0		Right Coronary
C1	5	65	0/038	0		Abdominal Vessels
C2	5	65	0/038	0		
SIM1	5	100	0/038	0		Intracranial vessels Abdominal Aortic vessels
SIM2	5	100	0/038	0		
VERT	5	100	0/038	0		Vertebral artery

PIG (S)	5	110	0/038	6		Aorta
PIG (S)	6	110	0/038	6		